

Title: Preparation of copper(II) chloride

Work instructions

Task: React 3 g of copper with hydrochloric and nitric acids to prepare $\text{CuCl}_2 \cdot 2\text{H}_2\text{O}$.

Theory

Copper(II) chloride is the second most abundant copper compound after CuSO_4 . Anhydrous is yellow brown, slowly absorbing moisture to form a green dihydrate.

It is commercially prepared by chlorination of copper:



It can also be prepared from CuO , $\text{Cu}(\text{OH})_2$ or $\text{Cu}(\text{CO}_3)_2$ by the action of hydrochloric acid.

CuCl_2 is used e.g. in organic syntheses or for pyrotechnical effects (it colours the flame blue green).

Equipment: graduated cylinder, glass rod, separating funnel, funnel, burner, ribbed funnel, Büchner funnel, suction flask, filter paper, watch glass, scales, porcelain bowl

Chemicals: hydrochloric acid (36%), nitric acid (63%), copper powder

Procedures:

1. To the weighed amount of Cu in the porcelain bowl, add the calculated amount (1.2 times the theory) of 24% HCl.
2. Cover the mixture with an inverted funnel and carefully add the calculated amount of 33% HNO_3 from the separating funnel, one part at a time, through the funnel stem.
3. When steam of NO_x stops forming, remove the funnel, rinse in a bowl, and heat the mixture gently until it dissolves.
4. Filter the solution obtained into a clean dish and evaporate with stirring. When a brown band of anhydrous CuCl_2 begins to form on the edge of the dish, wipe it off with a stick to prevent it from decomposing with heat. When the solution thickens, allow it to stand and aspirate the excluded crystals on a Büchner funnel.

Management of chemical substances

Chemicals	Form	H-statements	P-statements
HCl	Liquid, 36%	H290, H314, H335	P280, P303 + P361 + P353, P304 + P340, P305 + P351 + P338, P312
HNO ₃	Liquid, 65%	H290, H314	P260, P280, P303 + P361 + P353, P305 + P351 + P338, P310
Cu	Solid, powdery	H228, H315, H319, H335	P210, P273, P370 + P378
CuCl ₂ ·2H ₂ O	Solid	H290, H302 + H312, H315, H318	P302 + P352, P305 + P351 + P338, P321, P390, P501

Sources of risk and assessment of risk severity

Hydrochloric and nitric acids are strong inorganic acids whose fumes irritate the respiratory system. At the same time, reactions in which nitric acid acts as an oxidising agent often release NO_x. Therefore, work with them in a fume hood or a well-ventilated area at all times. Copper chloride is toxic to aquatic life, so avoid leaching it into sinks.

Waste management method

After weighing, transfer the product into the prepared container.

Risk reduction measures


Wear tight-fitting safety glasses or face shield, rubber gloves protective clothing and footwear. Wash hands thoroughly after handling. If skin contact occurs, wash with plenty of warm water and soap. In the event of an accident or if you feel unwell, inform the teacher immediately. Work in well-ventilated areas. Do not eat, drink, or smoke while working. Follow the safety instructions given by the teacher. If eyes are hit, rinse gently with water for a few minutes. Remove contact lenses, if fitted, and if they can be removed easily. Continue rinsing.

Name of the project: Digitization of chemistry experiments to improve the quality and support chemistry teaching in secondary schools
Acronym: ChemIQSoc
Project number: 2021-1-SK01-KA220-VET-000027995



Worksheet

Chemical calculations



Calculate:

- amount of 36% HCl + 1.2 times in g
- amount of 36% HCl in ml ($\rho_{\text{HCl},36\%} = 1.18 \text{ g/cm}^3$)
- quantity of H₂O for the preparation of 24% HCl solution
- amount of 63% HNO₃ + 1.2 times in g
- amount of 63% HNO₃ in ml ($\rho_{\text{HNO}_3,36\%} = 1.39 \text{ g/cm}^3$)
- quantity of H₂O for the preparation of 33% HNO₃ solution
- theoretical yield of CuCl₂·2H₂O
- relative yield of CuCl₂·2H₂O